



Recombinant Arabidopsis thaliana Uncharacterized mitochondrial protein AtMg00820 (AtMg00820)

Product Code	CSB-EP310097DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P92520
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MLTRSKAGIN KLNPKYSLTI TTTIKKEPKS VIFALKDPGW CQAMQEELDA LSRNKTWILV PPPVNQNILG CKWVFKTKLH SDGTLDRKA RLVAKGFHQE EGIYFVETYS PVRTATIRT ILNVAQQLEV GQSINWMFKM HFSGMIFKKK FICINLLVLR ILFIHPMCVC
Source	E.coli
Target Names	AtMg00820;
Protein Names	Recommended name: Uncharacterized mitochondrial protein AtMg00820 Alternative name(s): ORF170
Expression Region	1-170
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.