



Recombinant Arabidopsis thaliana Putative E3 ubiquitin-protein ligase SINAT1 (SINAT1)

Product Code	CSB-MP310126DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P93748
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MAPGGSALKE ALESNSTGVD YEVKMAKVEA NSKPTKSGSG SIGKFHSSNG VYELLECPVC TNLMPPIHQ CPNGHTLCSS CKLRVQNTCP TCRYELGNIR CLALEKVAES LEVPCRYQNL GCQDIFPYYS KLKHEQHCRF RSYSCP YAGS ECSVTGDIPT LVDHLKDDHK VDMHDGCTFN HRYVKS NPHE VENATWMLTV FNCFGRQFCL HFEAFQLGMA PVYMAFLRFM GDENEAKKFS YSLEVG AHSR KLTWQGIPRS IRDSHRKVRD SQDGLIIPRN LALYFSGSDK EELKLRVTGR IWKEE
Source	Mammalian cell
Target Names	SINAT1
Protein Names	Recommended name: Putative E3 ubiquitin-protein ligase SINAT1 EC= 6.3.2.- Alternative name(s): Seven in absentia homolog 1
Expression Region	1-305
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.