



Recombinant Solanum tuberosum 14-3-3-like protein 16R

Product Code	CSB-MP310127FIG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P93784
Product Type	Recombinant Protein
Immunogen Species	Solanum tuberosum (Potato)
Purity	>85% (SDS-PAGE)
Sequence	MASPREENVY MAKLAEQAER YEEMVEFMEK VVAAADGAEE LTVEERNLLS VAYKNVIGAR RASWRIISSI EQKEESRGNE DHVASIKEYR SKIESELTSI CNGILKLLDS KLIGSAATGD SKVFYLMKMG DYHRYLAEFK TGAERKEAAE NTLSAYKAAQ DIANAELAPT HPIRLGLALN FSVFYEILN SPDRACNLAK QAFDEAIAEL DTLGEESYKD STLIMQLLRD NLTLWTSDMQ DDGTDEIKEA APKPDNNE
Source	Mammalian cell
Protein Names	Recommended name: 14-3-3-like protein 16R
Expression Region	1-258
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.