



# Recombinant Phenylalanine dehydrogenase (pdh)

<b>Product Code</b>	CSB-YP310246SKM
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P97014
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Sporosarcina ureae
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MILVTLEQTL QDDKASVLDK MVEHEQILFC HDKATGLQAI IAVHDTTMGP ALGGCRMAY KTMDLALKDV LRLSKGMTYK CAAADVDFGG GKSVIIGDPL KDKTPEKFRA FGQFIESLNG RFYTGTDGMT TLEDFVHAMK ETNYIVGKPV EYGGGGDSSI PTALGVFYGI KATNQNLFGD DKVEGRKYSI QGLGKVGYKV AEHIINEGGN VIVTDINEQA IADIQKLGGS AVRVSSEEI YSQQADVFP CAFGGVINDD TLKVLKVRGI SGSANNQLAE SRHGELLREK GILYAPDYIV NGGGLIQVAD ELYGTNPARV LAKTENIYTS LLEV FHQAEQ DHMTTATAAD RMCEKRIADA KNRNSFFTQS NRPKWNFHQ
<b>Source</b>	Yeast
<b>Target Names</b>	pdh
<b>Protein Names</b>	Recommended name: Phenylalanine dehydrogenase Short name= PheDH EC= 1.4.1.20
<b>Expression Region</b>	1-379
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.