



# Recombinant Pig Pancreatic triacylglycerol lipase (PNLIP)

<b>Product Code</b>	CSB-EP018262PI-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P00591
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Sus scrofa (Pig)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	SEVCFPRLGC FSDDAPWAGI VQRPLKILPW SPKDVDRFL LYTNQNNY QELVADPSTI TNSNFRMDRK TRFIIHGFIK KGEEDWLSNI CKNLFKVESV NCICVDWKGK SRTGYTQASQ NIRIVGAEVA YFVEVLKSSL GYSPSNVHVI GHSLGSHAAG EAGRRTNGTI ERITGLDPAE PCFQGTPELV RLDPSDAKFV DVIHTDAAPI IPNLGFGMSQ TVGHLDFFPN GGKQMPGCQK NILSQVIDID GIWEGTRDFV ACNHLRSYKY YADSILNPDG FAGFPCDSYN VFTANKCFPC PSEGCPQMGH YADRFPGKTN GVSQVFLNT GDASNFRWR YKVSVTLSGK KVTGHILVSL FGNEGNSRQY EIKGTLQPD NTHSDEFDSD VEVGDLQVKV FIWYNNNVIN PTLPRVGASK IVERNDGKV YDFCSQETVR EEVLLTLNPC
<b>Source</b>	E.coli
<b>Target Names</b>	PNLIP
<b>Protein Names</b>	Recommended name: Pancreatic triacylglycerol lipase Short name= PL Short name= PTL Short name= Pancreatic lipase EC= 3.1.1.3
<b>Expression Region</b>	1-450
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full length protein
<b>Target Details</b>	This gene is a member of the lipase gene family. It encodes a carboxyl esterase that hydrolyzes insoluble, emulsified triglycerides, and is essential for the efficient digestion of dietary fats. This gene is expressed specifically in the pancreas.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.