



# Recombinant Pig Fructose-1,6-bisphosphatase 1 (FBP1)

<b>Product Code</b>	CSB-EP008459PI
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P00636
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Sus scrofa (Pig)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	TDQAAFDTN IVTLTRFVME EGRKARGTGE MTQLLNSLCT AVKAISTAVR KAGIAHLYGI AGSTNVTGDQ VKKLDVLSND LVINVLKSSF ATCVLVSEED KNAIIVEPEK RGKYVVC FDP LDGSSNIDCL VSIGTIFGIY RKNSTDEPSE KDALQPGRNL VAAGYALYGS ATMLVLAMVN GVNCFMLDPA IGEFILVDRD VKIKKKGSIY SINEGYAKEF DPAITEYIQR KKFPPDNSAP YGARYVGS MV ADVHRTL VYG GIFMYPANKK SPKGLRLLY ECNPMAYVME KAGGLATTGK EAVLDIVPTD IHQRAPIILG SPEDVTELLE IYQKHA AK
<b>Source</b>	E.coli
<b>Target Names</b>	FBP1
<b>Protein Names</b>	Recommended name: Fructose-1,6-bisphosphatase 1 Short name= FB Pase 1 EC= 3.1.3.11 Alternative name(s): D-fructose-1,6-bisphosphate 1- phosphohydrolase 1
<b>Expression Region</b>	2-338
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	Fructose-1,6-bisphosphatase 1, a gluconeogenesis regulatory enzyme, catalyzes the hydrolysis of fructose 1,6-bisphosphate to fructose 6-phosphate and inorganic phosphate. Fructose-1,6-diphosphatase deficiency is associated with hypoglycemia and metabolic acidosis.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.