



Recombinant Putative GTP-binding protein tag-210 (tag-210)

Product Code	CSB-MP310489CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P91917
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MPPKKNVVEE KPALIGRLGT NLKVGILGLP NVGKSTFFNV LTKSEAQAEN FPFCTIDPNE SRVAVQDDRF DWLVNHYKPA SKVPAFLNVT DIAGLVKGAS EGQGLGNAFL SHVSACDALF HLCRAFDDDD VTHVEGEVDP VRDLEISNE LFAKDLQFID GPLDKVEKLF TRANDKTKKI EYDTLVRVKK CLEEKKPVRQ ELWNEKEIEI LNKHLFLTAK PIVYLVNLSE KDYIRKKNKW LPKIKAWIDT NDAGAVLIPF SGAFELKLLD MPEDERQKYL KEQGVTSNLD KIVHTGYKAL QLEYFFTSGE DEVKAWTIQV GTPAPKAAGR IHTDFEKGFI MAEVMKVADL IELGDEAKCK AGGKYRQQGK TYIVQDGDVI FFKFNAGAGL QAKKK
Source	Mammalian cell
Target Names	ola-1
Protein Names	Recommended name: Putative GTP-binding protein tag-210
Expression Region	1-395
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.