



Recombinant Arabidopsis thaliana 50S ribosomal protein L24, chloroplastic (RPL24)

Product Code	CSB-MP310506DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P92959
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	KLKRWERKEC KPNSLPILHK MHVKFGDTVK VISGRDKGKI GEVTKIFTHN STIVIKDVNL KTKHMK SREE GEPGQIVKIE APIHSSNVML YSKEKDVVSR VGHKVLEDGQ KVRYLKLTGE LIDTIEKWKL LKEAKDKETT QVAVTSAS
Source	Mammalian cell
Target Names	RPL24
Protein Names	Recommended name: 50S ribosomal protein L24, chloroplastic Alternative name(s): CL24
Expression Region	51-198
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.