



Recombinant *Oryza sativa* subsp. japonica Gibberellin 20 oxidase 1 (20ox1)

Product Code	CSB-MP310529OFG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P93771
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MSMVVQQEQE VVFDAAVLSG QTEIPSQFIW PAEESPGSVA VEELEVALID VGAGAERSSV VRQVGEACER HGFFLVVNHG IEAALLEEAH RCMDAFFTLP LGEKQRAQRR AGESCGYASS FTGRFASKLP WKETLSFRYS SAGDEEGEEG VGEYLVRKLG AEHGRRLEGEV YSRICHEMSR LSLELMEVLG ESLGIVGDRR HYFRRFFQRN DSIMRLNYYP ACQRPLDTLG TGPHCDPTSL TILHQDHVGG LEVWAEGRWR AIRPRPGALV VNVGDTFMAL SNARYRSCLH RAVVNSTAPR RSLAFFLCPE MDTVVRPPEE LVDDHHPRVY PDFTWRALLD FTQRHYRADM RLFQAFSDWL NHHRHLQPTI YS
Source	Mammalian cell
Target Names	20ox1
Protein Names	Recommended name: Gibberellin 20 oxidase 1 EC= 1.14.11.- Alternative name(s): GA 20-oxidase 1 Gibberellin C-20 oxidase 1 Os20ox
Expression Region	1-372
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.