



Recombinant Bovine Vitamin K-dependent protein C (PROC)

Product Code	CSB-EP018742BO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P00745
Product Type	Recombinant Protein
Immunogen Species	Bos taurus (Bovine)
Purity	>85% (SDS-PAGE)
Sequence	A NSFLEELRPG NVERECSEEV CEFEEAREIF QNTEDTMAFW SFYSDGDQCE DRPSGSPCDL PCCGRGKCID GLGGFRCDCA EGWEGRFCLH EVRFSNCSAE NGGCAHYCME EEGRRHCSCA PGYRLEDDHQ LCVSKVTFPC GRLGKRMEKK RKTLKRDTNQ VDQKDQLDPR IVDGQEAGWG ESPWQAVLLD SKKKLVCGAV LIHVSWVLTV AHCLDSRKKL IVRLGEYDMR RWESWEVDLD IKEVIIHPNY TKSTSDNDIA LLRLAKPATL SQTIVPICLP DSGLSERKLT QVGQETVVTG WGYRDETKRN RTFVLSFIKV PVVYPYNACVH AMENKISENM LCAGILGDPR DACEGDSGGP MVTFFRGTW F LVGLVSWGEG CGRLYNYGVY TKVSRYL DWI YGHIKAQEAP LESQVP
Source	E.coli
Target Names	PROC
Protein Names	Recommended name: Vitamin K-dependent protein C EC= 3.4.21.69 Alternative name(s): Anticoagulant protein C Autoprothrombin IIA Blood coagulation factor XIV Cleaved into the following 3 chains: 1. Vitamin K-dependent protein C li
Expression Region	40-456
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.