



Recombinant Mycobacterium tuberculosis 50S ribosomal protein L2 (rpIB)

Product Code	CSB-EP310574MVZ-B
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P9WHA5
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Mycobacterium tuberculosis (strain ATCC 25618 / H37Rv)
Purity	>85% (SDS-PAGE)
Sequence	MAIRKYKPTT PGRRGASVSD FAEITRSTPE KSLVRPLHGR GGRNAHGRIT TRHKGGGHGR AYRMIDFRRN DKDGVNAKVA HIEYDPNRTA RIALHLYLDG EKRYIIAPNG LSQGDVVESG ANADIKPGNN LPLRNIPAGT LIHAVELRPG GGAKLARSAG SSIQLLGKEA SYASLRMPSTG EIRRVDVRCR ATVGEVGNAE QANINWGKAG RMRWKGKRPS VRGVVMNPVD HPHGGGEGKT SGGRHPVSPW GKPEGRTRNA NKSSNKFIVR RRRTGKKHSR
Source	E.coli
Target Names	rpIB
Protein Names	Recommended name: 50S ribosomal protein L2
Expression Region	1-280
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.