



# Recombinant Probable phenylalanine-4-hydroxylase 1 (pah-1)

<b>Product Code</b>	CSB-EP310865CXY
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P90925
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Caenorhabditis elegans
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MPPAGQDDLD FLKYAMESYV ADVNADIGKT TIVFTLREKA GALAETLKLK QAHDVNLISHI ESRPSKTHEG CYEVLVEFAE AEDHRKIEGV IEHFQQKA EK KVLVQDWNTK NKQNKDSVPW FPQKINDIDQ FANRILSYGA ELDADHPGFK DMTYRERRKF FADIAFNFKH GDKIPTITYT DEEIATWRTV YNELTVMYPK NACQEFNYIF PLLQQNCGFG PDRIPQLQDV SDFLKDCTGY TIRPVAGLLS SRDFLAGLAF RVFHSTQYIR HHSAPKYTPE PDICHELLGH VPLFADVEFA QFSQEIGLAS LGAPDDVIEK LATLYWFTIE FGICQQDGEK KAYGAGLLSS FGELQYALSD KPEVVDFDPA VCCVTKYPT EYQPKYFLAE SFASAKNKLK SWAATINRPF QIRYNAYTQR VEILDKVAAL QRLARDIRSD ISTLEEALGK VNNLKMK
<b>Source</b>	E.coli
<b>Target Names</b>	pah-1
<b>Protein Names</b>	Recommended name: Probable phenylalanine-4-hydroxylase 1 Short name= PAH EC= 1.14.16.1 Alternative name(s): Biogenic amine synthesis protein 2 Phe-4-monooxygenase
<b>Expression Region</b>	1-457
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.