



Recombinant Probable phenylalanine-4-hydroxylase 1 (pah-1)

Product Code	CSB-MP310865CXY
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P90925
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MPPAGQDDLD FLKYAMESYV ADVNADIGKT TIVFTLREKA GALAETLKLK QAHDVNLSHI ESRPSKTHEG CYEVLVEFAE AEDHRKIEGV IEHFQQKAEK KVLVQDWNTK NKQNKDSVPW FPQKINDIDQ FANRILSYGA ELDADHPGFK DMTYRERRKF FADIAFNFKH GDKIPTITYT DEEIATWRTV YNELTVMYPK NACQEFNYIF PLLQQNCGFG PDRIPQLQDV SDFLKDCTGY TIRPVAGLLS SRDFLAGLAF RVFHSTQYIR HHSAPKYTPE PDICHELLGH VPLFADVEFA QFSQEIGLAS LGAPDDVIEK LATLYWFTIE FGICQQDGEK KAYGAGLLSS FGELQYALSD KPEVVDFDPA VCCVTKYPT EYQPKYFLAE SFASAKNKLK SWAATINRPF QIRYNAYTQR VEILDKVAAL QRLARDIRSD ISTLEEALGK VNNLKMK
Source	Mammalian cell
Target Names	pah-1
Protein Names	Recommended name: Probable phenylalanine-4-hydroxylase 1 Short name= PAH EC= 1.14.16.1 Alternative name(s): Biogenic amine synthesis protein 2 Phe-4-monooxygenase
Expression Region	1-457
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.