



# Recombinant Human Glutathione reductase, mitochondrial (GSR)

<b>Product Code</b>	CSB-BP009968HU
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P00390
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MACRQEP QPQGPPPAAG AVASYDYLVI GGGSGGLASA RRAAELGARA AVVESHKLGG TCNVVGCVPK KVMWNTAVHS EFMHDHADYG FPSCEGKFNW RVIKEKRDAY VSRLNAIYQN NLTKSHIEII RGHAAFTSDP KPTIEVSGKK YTAPHILIAT GGMPSTPHES QIPGASLGIT SDGFFQLEEL PGRSVIVGAG YIAVEMAGIL SALGSKTSLM IRHDKVLRFS DSMISTNCTE ELENAGVEVL KFSQVKEVKK TLSGLEVSMV TAVPGRLPVM TMIPDVDCLL WAIGRVPNTK DLNKLGIQ TDDKGHIIVD EFQNTNVKGI YAVGDVCGKA LLTPVAIAAG RKLAHRLFY KEDSKLDYNN IPTVVFSHPP IGTVGLTEDE AIHKYGIENV KTYSTSFTPM YHAVTKRRTK CVMKMCANK EEKVVGIHMQ GLGCDEMLQG FAVAVKMGAT KADFDNTVAI HPTSSEELVT LR
<b>Source</b>	Baculovirus
<b>Target Names</b>	GSR
<b>Protein Names</b>	Recommended name: Glutathione reductase, mitochondrial Short name= GR Short name= GRase EC= 1.8.1.7
<b>Expression Region</b>	44-522
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.