



Recombinant Arabidopsis thaliana

Uncharacterized mitochondrial protein AtMg01410 (AtMg01410)

Product Code	CSB-MP310894DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P92567
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	≥85% (SDS-PAGE)
Sequence	MFGGRRRRRLP QDGTFNQTQP FDRLVGSRHS FSFDLKSATD RWPLVFLFEV VQYLFDRSFA SSVVNSAFAC NIFEVPFVKL KRRFSQVCFV AGQPLGYHGS WPTFALSHHI LVWWCAKQVH PGVRFTSYAV LGDDVVIADQ EVAKVYESAL GGLGVKISYQ KSLIPIQVLL SLLNASGLGN LVLKREIFPR ESGVLSLMSV PSLA
Source	Mammalian cell
Target Names	AtMg01410
Protein Names	Recommended name: Uncharacterized mitochondrial protein AtMg01410 Alternative name(s): ORF204
Expression Region	1-204
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.