



Recombinant Arabidopsis thaliana Ras-related protein RABB1c (RABB1C)

Product Code	CSB-MP310899DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P92963
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MSYAYLFKYI IIGDTGVGKS CLLLQFTDKR FQPVHDLTIG VEFGARMITI DNKPIKLQIW DTAGQESFRS ITRSYYRGAA GALLVYDITR RETFNHLASW LEDARQHANA NMTIMLIGNK CDLAHRRRAVS TEEGEQFAKE HGLIFMEASA KTAQNVVEAF IKTAATIYKK IQDGVFDVSN ESYGIKVGYG GIPGPSGGRD GSTSQGGGCC G
Source	Mammalian cell
Target Names	RABB1C
Protein Names	Recommended name: Ras-related protein RABB1c Short name= AtRABB1c Alternative name(s): Ras-related protein Rab2A Short name= AtRab2A
Expression Region	1-211
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.