



Recombinant Arabidopsis thaliana Ran-binding protein 1 homolog c (RANBP1C)

Product Code	CSB-BP310901DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P92985
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MASTEPEREN REDETEVNED EDTGAQVAPI VRLEEVAVTT GEEDEDAVLD LKSKMYRFDK EGNQWKERGA GTVKLLKHKE TGKVRLVMRQ SKTLKICANH LISSGMSVQE HSGNEKSCLW HATDFSDGEL KDELFCIRFA SIENCKTFME KFTEIAESQQ VGKESTQGDE AAGLIENLSV EENISEEKAK EAEEKEPAKE DKETKKEKVE EEKKTEAST
Source	Baculovirus
Target Names	RANBP1C
Protein Names	Recommended name: Ran-binding protein 1 homolog c
Expression Region	1-219
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.