



Recombinant Arabidopsis thaliana Germin-like protein subfamily 1 member 20 (GLP5A)

Product Code	CSB-EP310902DOA-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P92996
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	YDPSPLQD FCVAIDDLKG VFNNGRFCKD PKRVDAKDFE FSGLNVPGNT NNQVGSNVTT VNVDQIPGLN TMGISLVRID YAPHGQNPPH THPRGSEILV LVEGTLYVGF VSSNQDNNRL FAKVLHPGDV FVFPIGMIHF QVNVGKIPAV AFAGLSSQNA GVITIAITVF GSNPPIYPEL LARAFQLDAS VVKELQAKFG SI
Source	E.coli
Target Names	GLP5A
Protein Names	Recommended name: Germin-like protein subfamily 1 member 20 Alternative name(s): GLP2a copy 2 GLP2b Germin type 2 Short name= At-GERM2 Short name= AtGER2
Expression Region	23-222
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.