



Recombinant Arabidopsis thaliana

Uncharacterized mitochondrial protein AtMg00300 (AtMg00300)

Product Code	CSB-MP310913DOA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P93293
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	MRKGERMLAT SKEWERPRPV EASCSEGVLK VLKGCRTILK GNRHDSLYIL QGSVETGESN LAETAKDETR LWHSRLAHMS QRGMELLVKK GFLDSSKVSS LKFCEDCIYG KTHRVNFSTG QHTTKNPLDY VHSDLWGAPS VPLSF
Source	Mammalian cell
Target Names	AtMg00300;
Protein Names	Recommended name: Uncharacterized mitochondrial protein AtMg00300 Alternative name(s): ORF1451 ORF145a
Expression Region	1-145
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.