



Recombinant Nicotiana tabacum Actin-104

Product Code	CSB-MP310918NHE
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P93375
Product Type	Recombinant Protein
Immunogen Species	Nicotiana tabacum (Common tobacco)
Purity	>85% (SDS-PAGE)
Sequence	AGFAGDDAPR AVFPSIVGRP RHQGVMVGMG QKDAYVGDEA QSKRGILTLK YPIEHGIVSN WDDMEKIWHH TFYNELRVAP EEHPVLLTEA PLNPKANREK MTQIMFETFN VPAMYVAIQA VLSLYASGRT TGIVLDSGDG VSHTVPIYEG YALPHAILRL DLAGRGLTDY LMKILTERGY MFTTSAEREI VRDVKEKLAY VALDFEQELE TAKSSSSVEK SYELPDGQVI TIGAERFRCP EVLFQPSFIG MEAAGIHETT YNSIMKSDVD IRKDLYGNIV LSGGSTMFPG IADRMSKEIT ALAPSSMKIK VVAPPERKYS VWIGGSILAS LSTFQQ
Source	Mammalian cell
Protein Names	Recommended name: Actin-104
Expression Region	1-336
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.