



Recombinant *Oryza sativa* subsp. japonica S-adenosylmethionine synthase 2 (SAM2)

Product Code	CSB-EP310920OFG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P93438
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	MAAETFLFTS ESVNEGHPDK LCDQVSDAVL DACLAQDPDS KVACETCTKT NMVMVFGEIT TKATVDYEKI VRDTCRGIGF VSDDVGLDAD RCKVLVNIQ QSPDIAQGVH GHFTKRPEEI GAGDQGHMFG YATDEPELM PLSHVLATKL GARLTEVRKN GTCAWLRPDG KTQVTVEYLN DAGAMVPVRV HTVLISTQHD ETVTNDIEIAA DLKEHVIKPV IPDKYLDEKT IFHLNPSGRF VIGGPHGDAG LTGRKIIIDT YGGWGAHGGG AFSGKDPTKV DRSGAYIARQ AAKSIVASGL ARRCIVQVSY AIGVPEPLSV FVDSYGTGKI PDKEILKIVK ENDFRPGMM TINLDLKRGG NRFIKTAAYG HFGREDPDFT WEVVKPLKYE KASS
Source	E.coli
Target Names	SAM2
Protein Names	Recommended name: S-adenosylmethionine synthase 2 Short name= AdoMet synthase 2 EC= 2.5.1.6 Alternative name(s): Methionine adenosyltransferase 2 Short name= MAT 2
Expression Region	1-394
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.