



Recombinant Solanum tuberosum Guanine nucleotide-binding protein subunit beta (GB1)

Product Code	CSB-MP310922FIG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P93563
Product Type	Recombinant Protein
Immunogen Species	Solanum tuberosum (Potato)
Purity	>85% (SDS-PAGE)
Sequence	MSVAELKERH MAATQTVNDL REKLKQKRLQ LLDTDVSGYA KRQ GKSPVTF GPTDLVCCRI LQGHTGKVYS LDWTPEKNRI VSASQDGR LI VWNALTSQKT HAIKLPCAWV MTCAFSPSGQ SVACGGLDSA CSIFNLNSPI DKDGIHPVSR MLSGHKGYVS SCQYVPDEDT HLITSSGDQT CVLWDITTGL RTSVFGGEFQ SGHTADVLSV SISSSNPKLF VSGSCDTTAR LWDTRVASRA QRTFHGHESD VNTVKFFPDG NFRGTGSDDG SCRLFDIRTG HQLQVYNQPH GDGDIPHVTS MAFSISGRLL FVGYSNGDCY VWDTLLAKVV LNLGSVQNSH EGRISCLGLS ADGSALCTGS WDTNLKIWAF GGHR SVV
Source	Mammalian cell
Target Names	GB1
Protein Names	Recommended name: Guanine nucleotide-binding protein subunit beta
Expression Region	1-377
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.