



Recombinant Mouse Proto-oncogene serine/threonine-protein kinase mos (Mos)

Product Code	CSB-BP014733MO
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P00536
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	>85% (SDS-PAGE)
Sequence	MWLVLRIKEE GKGTGEGSNL QPCSQTS LAV PTHFSLVSHV TVPSEGVMP PLSLCRYLPR ELSPVDSRS CSIPLVAPRK AGKFLGTTP PRAPGLPRRL AWFSIDWEQV CLMHRLGSGG FGSVYKATYH GVPVAIKQVN KCTKDLRASQ RSFWAELNIA RLRHDNIVRV VAASTRTPED SNSLGTIIME FGGNVTLHQV YIGATRSPEP LSCREQLSLG KCLKYSLDVV NGLLFLHSQS ILHLDLKPAN ILISEQDVCK ISDFGCSQKL QDLRCRQASP HHIGGTYTHQ APEILKGEIA TPKADIYSFG ITLWQMTTRE VPYSGEPQYV QYAVVAYNLR PSLAGAVFTA SLTGKTLQNI IQSCWEARAL QRPGAELLQR DLKAFRGALG
Source	Baculovirus
Target Names	Mos
Protein Names	Recommended name: Proto-oncogene serine/threonine-protein kinase mos EC= 2.7.11.1 Alternative name(s): Oocyte maturation factor mos Proto-oncogene c-Mos
Expression Region	1-390
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.