



Recombinant *Sulfolobus solfataricus* 30S ribosomal protein S4P (rps4p)

Product Code	CSB-EP310998FPM
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P95987
Product Type	Recombinant Protein
Immunogen Species	<i>Sulfolobus solfataricus</i> (strain ATCC 35092 / DSM 1617 / JCM 11322 / P2)
Purity	≥85% (SDS-PAGE)
Sequence	MGDPKKSRRK WETPGHPWIK ERIGYEQELL GKYGLRNKRE IWIAQSIIRK FRHQARSLLA LPPAERAVRE KQLVGKLLKM GLLKKTATV DDILSLTEQD LLERRLQTIY YKKGLSNTIY QARQLITHGH IAVNGKRVTS PGYIVNVDEE NLIDYYVTSS FKSRRPPVMSQ QEGGEIGVKQ A
Source	E.coli
Target Names	rps4
Protein Names	Recommended name: 30S ribosomal protein S4P
Expression Region	1-181
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.