



Recombinant Goat Alpha-lactalbumin (LALBA)

Product Code	CSB-EP012724GO-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P00712
Product Type	Recombinant Protein
Immunogen Species	Capra hircus (Goat)
Purity	>85% (SDS-PAGE)
Sequence	E QLTKCEVFQK LKDLKDYGGV SLPEWVCTAF HTSGYDTQAI VQNNNSTEYGF LQINNKIWC KDDQNPHSRN ICNISCDFKL DDDLTDDIVC AKKILDKVGI NYWLAHKALC SEKLDQWLCE KL
Source	E.coli
Target Names	LALBA
Protein Names	Recommended name: Alpha-lactalbumin Alternative name(s): Lactose synthase B protein
Expression Region	20-142
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes alpha-lactalbumin, a principal protein of milk. Alpha-lactalbumin forms the regulatory subunit of the lactose synthase (LS) heterodimer and beta 1,4-galactosyltransferase (beta4Gal-T1) forms the catalytic component. Together, these proteins enable LS to produce lactose by transferring galactose moieties to glucose. As a monomer, alpha-lactalbumin strongly binds calcium and zinc ions and may possess bactericidal or antitumor activity. A folding variant of alpha-lactalbumin, called HAMLET, likely induces apoptosis in tumor and immature cells.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.