



Recombinant Candida boidinii Peroxisomal membrane protein PMP30B (PEX11B)

Product Code	CSB-MP311087CZG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q00317
Product Type	Recombinant Protein
Immunogen Species	Candida boidinii (Yeast)
Purity	>85% (SDS-PAGE)
Sequence	MVYGELIYHP VVTKLLKFLD SSASREKLLR LLQYLCRFLT FYTFKRNFNI ETIQLIKKIQ SSIGISRKPL RFLKNLPHLK NLNKIYSNEL LDSTLKIGDL IKNFGYALYF QFDTLQWLKL LGLLTSKNSG SLYFKIDKLA ANFWLIGLTG SIITDLRNLK ISYDSNKALL NEINSQNNNS NNDTLDEKLI EQNNDLILKN NEKINLNKRD LFKNILD SLI ALKGSQ LIDL NDGVLGFAGI ITSII GIEDI WNATKA
Source	Mammalian cell
Target Names	PEX11B
Protein Names	Recommended name: Peroxisomal membrane protein PMP30B Alternative name(s): Peroxin-11B Peroxisomal membrane protein PMP32
Expression Region	1-256
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.