



# Recombinant Schizosaccharomyces pombe G1/S-specific cyclin pas1 (pas1)

<b>Product Code</b>	CSB-BP311224SXV
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P87049
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSKVNMRALL ESFSTVILSI YPLAKSTNQQ SHSHLSLSTQ FLYETIRRSR VDYTTLLLLAL YYFIRFRDAA SSKPTNYIPL LCGRRLFLVC LMAATKFLQD RSFSNRAWRSR LSGLPVDKLL VLEYMFYQCI DYRLVVPKHI FARWSLLVGE CCAHATARYD STDPNVASFG SVAQYWVSLF SNVQSSLLDDL FSIACLTkia HRRMNANAAL KNQATRKpSS SPQTTQDSSP ILTMAPSTPV SVGSTPPSTP SVLPIAKQLA PMNVCKAHIQ ASNQSRtlTT ASPPEQIPLM EPQVYVNPQV LPGRLSSLSK PVSLPPTPSS PKVGVYRPMT SKSNGGVAYC YNAQKLNNAT GPVTFNMPFA SVLPLAVSVS CDLGTASAYV ASLPQPCSQK RHLEEDYSCL TEHSAKRrsy F
<b>Source</b>	Baculovirus
<b>Target Names</b>	pas1
<b>Protein Names</b>	Recommended name: G1/S-specific cyclin pas1
<b>Expression Region</b>	1-411
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.