



Recombinant Schizosaccharomyces pombe Uba3-binding protein but2 (but2)

Product Code	CSB-EP311232SXV-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P87167
Product Type	Recombinant Protein
Immunogen Species	Schizosaccharomyces pombe (strain 972 / ATCC 24843) (Fission yeast)
Purity	≥85% (SDS-PAGE)
Sequence	<pre>MQLLNSFLGF AASIAVLASS ADAAPTLYKR KSKSNTNTAK SVAILDGVGV NNTFGVISYR EKFPMDYHIW HTAPSTGKTY LQPWNDAVEP STYYIDEGTF LRTGLKFAYI GDNMDLAFTN HTDWKASKHD DGTHRVDLSQ RKPANNFMAT QRCGKLD PFG YQLKRSKKVF QACGTQVFVG SGRSIDCQWM DSVALNLSRY YGNTEALSSL PLPRNLSADI PAKSSRLFPH GIYNFNFSAP NKRMQRFTAS EATTVNGQNQ SVYLTYDVPY GRSLYYYTSC ALEFRFTNEF FPMDVTPNDG SAQFVVYNMS GNPKKQTTSS KGPTRLYEVA RFNCTTRGCE YTQNIPCRA GHSHTYELAP ASPNTSISWI QSYSPRIGVT LQVYSNANFD</pre>
Source	E.coli
Target Names	but2
Protein Names	Recommended name: Uba3-binding protein but2
Expression Region	1-390
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.