



Recombinant Drosophila melanogaster Protein stand still (stil)

Product Code	CSB-MP311271DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P92189
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	≥85% (SDS-PAGE)
Sequence	MSVNLDGSEE KMELKKHVPY ILNGELYRIE KQVGDNVTVK CCYCPPDRIY RGSVRSTGNF HMHIKRRHSS LLGKLHEMKV AALEERRDRI MKNRRFAKSR KKAPVAVAAT STAAQSDSGV FVDMQAAVPS GNESHELKIK TVFQRHKQEQ EGATRKLEDS TSDKANLPNI PKNLGIVVQN VSNISVETLP GGSTPASVSF LGRPVKPEQG SGPSFLADQP AAIDLSQVPP VQGESKSSGS LASSMEDVSM EYSRSQALSQ SL SMAHFLEH PQRDVLQRLE RSMAQISQEL HCRNRIEHNR MLLEAAKFKF LNPNFQFEPN L
Source	Mammalian cell
Target Names	stil
Protein Names	Recommended name: Protein stand still
Expression Region	1-321
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.