



Recombinant Mouse Four and a half LIM domains protein 1 (Fhl1)

| | |
|--------------------------|---|
| Product Code | CSB-EP008663MO-B |
| Storage | Store at -20°C, for extended storage, conserve at -20°C or -80°C. |
| Uniprot No. | P97447 |
| Product Type | Recombinant Protein |
| Immunogen Species | Mus musculus (Mouse) |
| Purity | >85% (SDS-PAGE) |
| Sequence | SEKFDCHYC RDPLQGKKYV QKDGRHCCLK CFDKFCANTC VDCRKPI SAD AKEVHYKNRY WHDNCFRCAK CLHPLASETF VSKDGGKILCN KCATREDSPR CKGCFKAIVA GDQNV EYKGT VWHKDCFTCS NCKQVIGTGS FFPKGEDFYC VTCHETKFAK HCVKCNKAIT SGGITYQDQP WHAECFVCVT CSKLAGQRF TAVEDQYYCV DCYKNFVAKK CAGCKNPITG FGKGSSVVAY EGQSWHDYCF HCKKCSVNLA NKRFV FHNEQ VYCPDCAKKL |
| Source | E.coli |
| Target Names | Fhl1 |
| Protein Names | Recommended name: Four and a half LIM domains protein 1 Short name= FHL-1 Alternative name(s): KyoT RBP-associated molecule 14-1 Short name= RAM14-1 Skeletal muscle LIM-protein 1 Short name= SLIM Short name= SLI |
| Expression Region | 2-280 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | Full Length of Mature Protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |