



Recombinant *Saccharomyces cerevisiae* Cold sensitive U2 snRNA suppressor 1 (CUS1)

Product Code	CSB-YP311614SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q02554
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MARTKSRKRS GNNQNKNASV VNNKAEIAAM IDARRLEQKK KGGVTNSK GK TNKVVD AKLE KEFKDVLQRF QVQENDTPKE ITKDEKNNHV VIVEKNPVMN RKHTAEDELE DTPSDGIEEH LSARKRRKTE KPSLSQLKSQ VPYPQII EWY DCDARYPGLL ASIKCTKNVI PVP SHWQSKK EYLSGRSLLG KRPFELPDII KKTNIEQMRS TLPQSGLDGQ DEKSLKEASR ARVQPKMGAL DLDYKKLHDV FFKIGANWKP DHLLCFGDVY YENRNLFEET NWKRMVDH KR PGRISQELRA IMNLPEGQLP PWCMKMKDIG LPTGY PDLKI AGLNWDITNL KGDVYGKIIP NHHSRSKKQG RNYFGALISF ETPEFENSKE DTQANAENGR QDDKIDDEVE HKLDHFQEDI SEVTS AEEKL ERNEEESEKQ LYTVLK
Source	Yeast
Target Names	CUS1
Protein Names	Recommended name: Cold sensitive U2 snRNA suppressor 1
Expression Region	1-436
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.