



Recombinant *Saccharomyces cerevisiae*

Uncharacterized protein YMR265C (YMR265C)

Product Code	CSB-BP311792SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q03508
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MEEFEEFRRK GEMSSRCGNH RVLRKWNSCA CELAVPFEVP EHAITKLHIY DFDNTLFATP GPTEQLYTRE LLNLLTSSTL PNGGWWNEPG FLQAAIEISK TKPRRYSWNA DIVKLAEESEY SAKDTISIVL TGREESKFHK LIEHALQTAR SHWKCSENEFRFNAVCLKKR AISEYTSKYK KELMRDFLEY YPSLRELSIY DDRIHQIDAF KSFFHSLDLP RLKWSAIPVR PFTKSLPREQ ELEMVMDMVR KNNSQALSTS QKFDLRRTPR QIGYILCTAS HRLLSIEVIK YLKRRKGRRT FRPKLYEHPL YIPCAEPGKD IPALEIAKWV SNNDTRTFDS EKKVQHISQI FYLEQPGKCI VHFQVTDLAV IASAHNRRK PLEVYFKATP EPNRYTFTLF PEYIVTGHFY KRDRIEDLEV VTERLINCKE DIHWVPLDNT IPIKAFFGRF AKLAAIPCSN A
Source	Baculovirus
Target Names	YMR265C
Protein Names	Recommended name: Uncharacterized protein YMR265C
Expression Region	1-461
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.