



Recombinant Drosophila melanogaster Cysteine string protein (Csp)

Product Code	CSB-BP311815DLU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q03751
Product Type	Recombinant Protein
Immunogen Species	Drosophila melanogaster (Fruit fly)
Purity	>85% (SDS-PAGE)
Sequence	MSAPGMDKRRK LSTSGDSLIE ILGLPKTATG DDIKKTYRKL ALKYHPDKNP DNVDAADKFK EVNRAHSILS NQTKRNIYDN YGSLGLYIAE QFGEENVNAY FVVTSPAVKA VVICCAVITG CCCCCCCCCC CNFCCGKFKP PVNESHQYS HLNRPDGNRE GNDMPHTLGQ PPRLEDVDLD DVNLGAGGAP VTSQPREQAG GQPVFAMPPP SGAVGVNPFT GAPVAANENT SLNTTEQTTY TPDMVNQKY
Source	Baculovirus
Target Names	Csp
Protein Names	Recommended name: Cysteine string protein
Expression Region	1-249
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.