



Recombinant *Saccharomyces cerevisiae* Trafficking protein particle complex subunit 23 (TRS23)

Product Code	CSB-EP311817SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q03784
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MAIETILVIN KSGGLIYQRN FTNDEQKLNS NEYLILASTL HGVFAIASQL TPKALQLTQQ TNIENTIPYI PYVGMSSNRS DTRNGGGNNN KHTNNEKLG FKGDDFFKEP FTNWNKSGLR QLCTDQFTMF IQTLTGLKF VAISSSVMPQ RQPTIATTDK PDRPKSTSNL AIQIADNFLR KVYCLYSDYV MKDPSYSMEM PIRSNLFDEK VKKMVENLQ
Source	E.coli
Target Names	TRS23
Protein Names	Recommended name: Trafficking protein particle complex subunit 23 Short name= TRAPP subunit 23 Alternative name(s): Transport protein particle 23 kDa subunit
Expression Region	1-219
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.