



Recombinant *Saccharomyces cerevisiae* Uncharacterized protein YMR111C (YMR111C)

Product Code	CSB-BP311995SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q04461
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MPAREYNYVE GFGGYGSLLD DDSDRDSERR NHDLGQRTIT TSPTGVS RHA ALNRYMIPGR INPLFRPTDA AQPIIVSTST SASATEPTNR IGPGRKETP ETNFNAFLIA QLTRMEEQNA NLKEEISLMK KEQELFFLEN QKKLEKGFKD INKYVEDVSA MKEVFKEVVG IMTGERIRFI DHTGENVTPQ EAARVGNPST STQAHQSQSR STNWQEYSMH ASILAGDPRI KPEPGLSDFE NGEYDGNESD ENATTRNLPL NNPDSVSNAD DSNQLDGTG NENDIRNRRG CVGTSYKLN R AIQNVTDAAAR EYFEGLPGQP SVLSLERRYG STWRRSAKER TLFTKRMTII KRIIDIKDDP SKYGLSLPEN KISRNQAIKV VENIRLGNNT FKGHHCLSM SQLYEYFSKK MDKLEDYSLT LKRRGKPRRI FLEEREARL SLQQPHSIPN SSTGTPEHDQ DT
Source	Baculovirus
Target Names	YMR111C
Protein Names	Recommended name: Uncharacterized protein YMR111C
Expression Region	1-462
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.