



# Recombinant Chicken Troponin C, skeletal muscle (TNNC2)

<b>Product Code</b>	CSB-YP024010CH
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P02588
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Gallus gallus (Chicken)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	ASMTDQQAE ARAFLSEEMI AEFKAAFDMF DADGGGDIST KELGTVMRML GQNPTKEELD AIIEEVEDG SGTIDFEEFL VMMVRQMKED AKGKSEEELA NCFRIFDKNA DGFIDIEELG EILRATGEHV TEEDIEDLMK DSDKNNDGRI DFDEFLKMME GVQ
<b>Source</b>	Yeast
<b>Target Names</b>	TNNC2
<b>Protein Names</b>	Recommended name: Troponin C, skeletal muscle
<b>Expression Region</b>	2-163
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	Troponin (Tn), a key protein complex in the regulation of striated muscle contraction, is composed of 3 subunits. The Tn-I subunit inhibits actomyosin ATPase, the Tn-T subunit binds tropomyosin and Tn-C, while the Tn-C subunit binds calcium and overcomes the inhibitory action of the troponin complex on actin filaments. This protein is the Tn-C subunit.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.