



Recombinant Macaca mulatta (Rhesus macaque) Apolipoprotein A-II (APOA2)

Product Code	CSB-BP001915MOW
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P02653
Product Type	Recombinant Protein
Immunogen Species	Macaca mulatta (Rhesus macaque)
Purity	>85% (SDS-PAGE)
Sequence	QAEPSVESL VSQYFQTVTD YGKDLMEKVK SPELQAQAKA YFEKSKEQLT PLVKKAGTDL VNFLSYFVEL RTQPATQ
Source	Baculovirus
Target Names	APOA2
Protein Names	Recommended name: Apolipoprotein A-II Short name= Apo-AII Short name= ApoA-II Alternative name(s): Apolipoprotein A2 Cleaved into the following chain: 1. Truncated apolipoprotein A-II
Expression Region	1-77
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full length protein
Target Details	This gene encodes apolipoprotein (apo-) A-II, which is the second most abundant protein of the high density lipoprotein particles. The protein is found in plasma as a monomer, homodimer, or heterodimer with apolipoprotein D. Defects in this gene may result in apolipoprotein A-II deficiency or hypercholesterolemia.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.