



Recombinant Rat Alpha-fetoprotein (Afp)

Product Code	CSB-BP001421RA
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P02773
Product Type	Recombinant Protein
Immunogen Species	Rattus norvegicus (Rat)
Purity	≥85% (SDS-PAGE)
Sequence	RVLHTN EFGIESTLDS SQCPTEKNMF NVATIVVAQF VQDATKAEVN KMSSDALAAM KENTGDGCLE NQLSVFLDEI CHETELSNKY GFSGCCNQSG VERHQCLLAR KKTAPDSVPP FHFETAESC PAYEENRAMS INTFIYDVSK RNPFLYAPTI LYLAAQYDKA VPACCKADNM EECFQTKRAS MAKELREGSM LNEHVCAVIR KFGSRNLQAV LIKLSQKFP KANITEIRKL ALDVAHIHEQ CCHGNAMECL QDGESVMTHM CSQQEILSSK TAECCKLPTI ELGYCIIHAE NGDKPEGLTL NPSEFLGDRN FAQFSSEEKL LFMASFLHEY SRNHPNLPVS VILKTAQSYQ EILEKCSQSE TPSKCQDNME EELQKHIQES QALAKQSCNL YQKLGPPYYLQ NLFLIGYTRK APQLTSAELI DLTGKMVSIA STCCQLSEEK RSACGEGLAD IYIGHLCLRH EANPVNSGIN HCCSSSYSNR RLCITSFLRD ETYVPPPFSE DKFIFHKDLC QAQGRALQTM KQELLINLVK QKPEMTEEQH AAVTADFSGL LEKCKDQDQ EACFAKEGPK LISKTREALG V
Source	Baculovirus
Target Names	Afp
Protein Names	Recommended name: Alpha-fetoprotein Alternative name(s): Alpha-1-fetoprotein Alpha-fetoglobulin
Expression Region	25-611
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes alpha-fetoprotein, a major plasma protein produced by the yolk sac and the liver during fetal life. Alpha-fetoprotein expression in adults is often associated with hepatoma or teratoma. However, hereditary persistence of alpha-fetoprotein may also be found in individuals with no obvious pathology. The protein is thought to be the fetal counterpart of serum albumin, and the alpha-fetoprotein and albumin genes are present in tandem in the same transcriptional orientation on chromosome 4. Alpha-fetoprotein is found in monomeric as well as dimeric and trimeric forms, and binds copper, nickel, fatty acids and bilirubin. The level of alpha-fetoprotein in amniotic fluid is used to measure renal loss of protein to screen for spina bifida and anencephaly.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a



concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.

Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.