



Recombinant *Saccharomyces cerevisiae* Chromosome stability protein 9 (CST9)

Product Code	CSB-MP312321SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q06032
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MGGYLAIVFI PQTNTKSMRE KKQKCLKQVR RLSLISPKKY IMPDSIFEQP FVYCGVCHRR TSHGDPLRLT SCAHILCSQH SPLTSKVCPI CRSSDISIIN LVESKQLPTD IRIFFEPLPP LLESLYNVSQ FQLNGLSKQC QYYQNHCLKL REKCARQQQL LYQAKIELDS MAILKKRIQE LESVLNHNHV SSMSVGVLP RNSHQNHYP PPTVDLTVDD NSLEEFKAKS FIKKLLKKNSS LRNSSKNNNG TVTPSTSGRV NKNQPLFMET LNNPNRNSIP PPGMNPANNS NLPNISTIAE STNLNRFSFS PVRVAKGFDG KLPNLDILT NGSVSSKNIS RLSSASLQPS SPLSSSSNRL ILPNSNLKEL HHSNTPLTST STQFPKSALEK LKITRKRNT ISGSNRITHN LSSHVRSSGL AFSSSSNSLQ QSKLPKSNIL KRSNSTQQLT NTHLKSDNHL PPRSSNTVLG SSKKNNKFRR IR
Source	Mammalian cell
Target Names	CST9
Protein Names	Recommended name: Chromosome stability protein 9 Alternative name(s): Molecular zipper protein 3
Expression Region	1-482
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.