



Recombinant *Saccharomyces cerevisiae* Protein RIF2 (RIF2)

Product Code	CSB-EP312340SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q06208
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MEHVDSDFAP IRRSKKVVDS DKIVKAISDD LEQKNFTVLR KLNLVPIKKS VSSPKVCKPS PVKERVDHVF YQKFKSMALQ ELGTNYLSIS YVPSLSKFLS KNLRSMKNCI VFFDKVEHIH QYAGIDRAVS ETLSLVDINV VIEMNDYLM KEGIQSSKSK ECIESMGQAS YSGQLDFEAS EKPSNHTSDL MMMVMRKINN DESIDHIVYF KFEQLDKLST STIIEPSKLT EFINVLSVLE KSNNIAFKVL IYSNNVSISS LLSTSLKKKL NTKYTVFEMP ILTCAQEQEY LKKMIKFTFD SGSKLLQSYN SLVTCQLNNK ESNLAIFFEF LKVFPHPFTY LFNAYTEIIV QSRTFDELLED KIRNRLTIKN YPHSAYNFKK NQRLPLKLTR KVHDR
Source	E.coli
Target Names	RIF2
Protein Names	Recommended name: Protein RIF2 Alternative name(s): RAP1-interacting factor 2
Expression Region	1-395
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.