



Recombinant *Saccharomyces cerevisiae* Uncharacterized protein YDR370C (YDR370C)

Product Code	CSB-YP312354SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q06349
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSTEQDAVLG LAKDLEGINL LTVPNLERGH QSKLCKEKTT SDSSSSRKPS QQRDNRYRKR PKLICIPYTS FLHTGMHNFL TKPPRDIFHE SKEVALFTNG RAYTILRKDL IPNLKESIAE LYESSLLEAK KRKVPYLGHD LFANIDEFVP MTISELDSVS PCFSYIENWI LDNPGKDFKI GKKFTVVTTR HHIVDLTMHL FNRRNRQTSL IVTYMGAGLL SFCRNVKKDS QMSKEGIYSN DPNMKKICYS GFEFENWVTE NSKVADLTGS KCPIFSLVES KLSEEIGLLI RCEMDAFNPV SETNTEKCF APLSMHNSNH RRKLLKTWVQ TGLLPNSDIM IGLRDSHSGQ LLDIQWYSRD LLCKKFNHPG LPTNKKELNY NAQIAVEWCH YCIEAICKLV EANISDYSST KPESFEIGID TNNAIVITKL KTTPRNVELF GM
Source	Yeast
Target Names	DXO1
Protein Names	Recommended name: Uncharacterized protein YDR370C
Expression Region	1-442
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.