



# Recombinant *Saccharomyces cerevisiae* Pre-mRNA-splicing factor URN1 (URN1)

<b>Product Code</b>	CSB-MP312373SVG
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	Q06525
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MRGEWQEFKT PAGKKYYYNK NTKQSRWEKP NLKKGSNLES NAKESQTERK PTFSLVLVNG WHLIIYNDGT KLYFNDDSK EFKNDISQEDD SRCRSLIESL DKEKLVLLIG VARGYTMREE DIDKILESCN EEIHLFKRNQ DEVERKDEIS EEAGDVKSPL QESHTGLVSG YGSSSGEED EEEDEEEDEN EEQIVNQDIS IIDDLNRIDT DDIDERNIFF ELFDYKLDK FSTWSLQSKK IENDPDFYKI RDDTVRESLF EEWCGERSGN ATAEESDSED NSEDDSEVLE PTKYHYLAQI VANAGTIAPD TIPQDIRKQQ KALYKAYKIK EYIPSKRDQD KFVSQLLFYI KTFDLEQRKE IFCDCLRDHE RDFTGAVESL RQDKELIDRW QTLLKAPADS SSIEDILLSI EHRCCVSPIV VTEPRYYVVG ILEKTVVWVR WLAAEVPSS RFTPVGAGNE PINPE
<b>Source</b>	Mammalian cell
<b>Target Names</b>	URN1
<b>Protein Names</b>	Recommended name: Pre-mRNA-splicing factor URN1 Alternative name(s): U2-U5-U6 snRNP, RES complex and NTC-interacting pre-mRNA-splicing factor 1
<b>Expression Region</b>	1-465
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.