



Recombinant *Saccharomyces cerevisiae* Uncharacterized protein YPR078C (YPR078C)

Product Code	CSB-MP312395SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q06813
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	<p>MQTISGVLPT VLSPSELRSD DERTFQFDEE AEITTHLTES EDLRRLINET AQLGVRVDHI HDKTDQEIAR LEKVIKEVTE SDTFFRSCSG WFKTNKNFSD SESSNTQLK SLSQLHGRYD RDWRQRLNKW FRKNKSKLAL PSDNNLEEVN DDKVYGYGED LMERGKTPYF SDIDDFMNGI NIISPLTPDD FENDDTLVKI DETCQIHSAS EPEKTSISPT FGKNIKKELV TDDTESIISG PPLQENKCTL LKYRYVRTSL DMLGSEKSSS KNNSSGGMFRI FHKSANFGDK NQENVPRVWD TLRNNLGREI YLLQGRFKKW TTKHQNLKKG QPCKDEDAVT VPLPSSDPGK ETQLETCLCF VPEPGDQPLV QA</p>
Source	Mammalian cell
Target Names	YPR078C
Protein Names	Recommended name: Uncharacterized protein YPR078C
Expression Region	1-372
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.