



Recombinant *Saccharomyces cerevisiae* Alcohol-sensitive RING finger protein 1 (ASR1)

Product Code	CSB-MP312397
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q06834
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	≥85% (SDS-PAGE)
Sequence	MEECPICLAD DQEGEQFGCL NVCGHKFHLN CIREWHKYSI NLKCPICRVE STHLEVGEQ HALSINLKMG FMIKNAIDYV GAETTNERNE DDTGEQDQEI EFLSERLRGT LVMDTIKIIQ CSICGDTDVS RLSLYCQDCE AIYHETCLRG LACEVGDRNT WQECTDCRSN ALLELRMGAI SSQLASYDSR NSMIFAGELR DKHSVKTQQM YEQIRNAKHK IQMHVRRALD RYPLPLLRFK DAYKHAVNKQV SRKLYRLSDN KYLPDQYDYD SLARTGVHTE LLIYCHDE
Source	Mammalian cell
Target Names	ASR1
Protein Names	Recommended name: Alcohol-sensitive RING finger http://www.uniprot.org/uniprot/Q00327
Expression Region	1-288
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.