



Recombinant Mouse Adenosine deaminase (Ada)

Product Code	CSB-EP001268MO-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P03958
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	AQTPAFNKP KVELHVHLDG AIKPETILYF GKKGIALPA DTVEELRNII GMDKPLSLPG FLAKFDYYMP VIAGCREAIK RIAYEFVEMK AKEGVVYVEV RYSPELLANS KVDPMPWNQT EGDVTPDDVV DLVNQGLQEG EQAFGIKVRS ILCCMRHQPS WSLEVLELCK KYNQKTVVAM DLAGDETI EG SSLFPGHVEA YEGAVKNGIH RTVHAGEVGS PEVVREAVDI LKTERVGHGY HTIEDEALYN RLLKENMHFE VCPWSSYLTG AWDPKTTHAV VRFKNDKANY SLNTDDPLIF KSTLTDYQM TKKDMGFTEE EFKRLNINAA KSSFLPEEEK KELLERLYRE YQ
Source	E.coli
Target Names	Ada
Protein Names	Recommended name: Adenosine deaminase EC= 3.5.4.4 Alternative name(s): Adenosine aminohydrolase
Expression Region	2-352
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	This gene encodes an enzyme that catalyzes the hydrolysis of adenosine to inosine. Various mutations have been described for this gene and have been linked to human diseases. Deficiency in this enzyme causes a form of severe combined immunodeficiency disease (SCID), in which there is dysfunction of both B and T lymphocytes with impaired cellular immunity and decreased production of immunoglobulins, whereas elevated levels of this enzyme have been associated with congenital hemolytic anemia.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.