



Recombinant Human Tissue alpha-L-fucosidase (FUCA1)

Product Code	CSB-EP009062HU-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P04066
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	<p>QPPRRYTPD WPSLDSRPLP AWFDEAKFGV FIHWGVFSVP AWGSEWFWWH WQGEGRPQYQ RFMRDNYPPG FSYADFGPQF TARFFHPEEW ADLQAAGAK YVVLTTKHHE GFTNWPSPVS WNWNSKDVGP HRDLVGELGT ALRKRNIYR YG LYHSLLEWFH PLYLLDKKNG FKTQHFVSAK TPELYDLVN SYKPDLIWSD GEWECPTYW NSTNFLSWLY NDSPVKDEVV VDRWGWQNC CHHGGYYNCE DKFKPQSLPD HKWEMCTSID KFSWGYRRDM ALSDVTEESE ISELVQTVS LGGNYLLNIG PTKDGLIVPI FQERLLAVGK WLSINGEAIY ASKPWRVQWE KNTTSVWYTS KGSVYAIFL HWPENGVNL ESPITTSTK ITMLGIQGD L KWSTDPDKGL FISLPQLPPS AVPAEFAWTI KLTGVK</p>
Source	E.coli
Target Names	FUCA1
Protein Names	Recommended name: Tissue alpha-L-fucosidase EC= 3.2.1.51 Alternative name(s): Alpha-L-fucosidase I Alpha-L-fucoside fucohydrolase 1 Short name= Alpha-L-fucosidase 1
Expression Region	32-466
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.