



Recombinant Human Prorelaxin H2 (RLN2)

Product Code	CSB-YP019750HU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	P04090
Storage Buffer	Lyophilized from Tris/PBS-based buffer, 6% Trehalose, pH 8.0
Product Type	Recombinant Proteins
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	DSWMEEVIKLCGRELVRAQIAICGMSTWSKRSLSQEDAPQTPRPVAEIVPSFIN KDTETINMMSEFVANLPQELKLTSEMQPALPQLQQHVPVLKDSSLLFEEFKKL IRNRQSEAADSSPSELKYLGLDTHSRKKRQLYSALANKCCHVGCTKRSLARFC
Source	Yeast
Target Names	RLN2
Protein Names	Recommended name: Prorelaxin H2 Cleaved into the following 2 chains: 1. Relaxin B chain 2. Relaxin A chain
Expression Region	25-185aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.