



Recombinant *Saccharomyces cerevisiae* U3 small nucleolar RNA-associated protein 6 (UTP6)

Product Code	CSB-MP312534SVG
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	Q02354
Product Type	Recombinant Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain ATCC 204508 / S288c) (Baker's yeast)
Purity	>85% (SDS-PAGE)
Sequence	MSKTRYYLEQ CIPEMDDLVE KGLFTKNEVS LIMKKRTDFE HRLNSRGSSI NDYIKYINYE SNVNKLRAKR CKRILQVKKT NSLSDWSIQQ RIGFIYQRGT NKFPQDLKFW AMYLNMYKAR GNQTSYKKIH NIYNQLLKLH PTNVDIWISC AKYEYEVHAN FKSCRNIFQN GLRFNPDVVK LWYEVVKFEL NFITKLINRR KVMGLINERE QELDMQNEQK NNQAPDEEKS HLQVPSTGDS MKDKLNELPE ADISVLGNAE TNPALRGDIA LTIFDVCMT LGKHYINKHK GYYAISDSKM NIELNKETLN YLFSESLRYI KLFDEFDLE RDYLINHVLQ FWKNDMYDLS LRKDLPELYL KTVMIDITLN IRYMPVEKLD IDQLQLSVKK YFAYISKLDS ASVKSLKNEY RSYLQDNYLK KMNAEDDPY KILDLIISKL
Source	Mammalian cell
Target Names	UTP6
Protein Names	Recommended name: U3 small nucleolar RNA-associated protein 6 Short name= U3 snoRNA-associated protein 6 Alternative name(s): U three protein 6
Expression Region	1-440
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.