



# Recombinant Mouse Fatty acid-binding protein, adipocyte (Fabp4)

<b>Product Code</b>	CSB-EP007945MO-B
<b>Storage</b>	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
<b>Uniprot No.</b>	P04117
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Mus musculus (Mouse)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	CDAFVGTWK LVSSNFDDY MKEVGVGFAT RKVAGMAKPN MIISVNGDLV TIRSESTFKN TEISFKLQVE FDEITADDRK VKSIITLDGG ALVQVQKWDG KSTTIKRKRD GDKLVVECVK KGVSTSTRVYE RA
<b>Source</b>	E.coli
<b>Target Names</b>	Fabp4
<b>Protein Names</b>	Recommended name: Fatty acid-binding protein, adipocyte Alternative name(s): 3T3-L1 lipid-binding protein Adipocyte lipid-binding protein Short name= ALBP Adipocyte-type fatty acid-binding protein Short name= A-FABP Short
<b>Expression Region</b>	2-132
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Target Details</b>	FABP4 encodes the fatty acid binding protein found in adipocytes. Fatty acid binding proteins are a family of small, highly conserved, cytoplasmic proteins that bind long-chain fatty acids and other hydrophobic ligands. It is thought that FABPs roles include fatty acid uptake, transport, and metabolism.
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.