



Recombinant Mouse Fatty acid-binding protein, adipocyte (Fabp4)

Product Code	CSB-EP007945MO-B
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P04117
Product Type	Recombinant Protein
Immunogen Species	Mus musculus (Mouse)
Purity	≥85% (SDS-PAGE)
Sequence	CDAFVGTWK LVSSSENFDDY MKEVGVGFAT RKVAGMAKPN MIISVNGDLV TIRSESTFKN TEISFKLGE FDEITADDRK VKSIITLDGG ALVQVQKWDG KSTTIKRKRD GDKLVVECV M KGV TSTRVYE RA
Source	E.coli
Target Names	Fabp4
Protein Names	Recommended name: Fatty acid-binding protein, adipocyte Alternative name(s): 3T3-L1 lipid-binding protein Adipocyte lipid-binding protein Short name= ALBP Adipocyte-type fatty acid-binding protein Short name= A-FABP Short
Expression Region	2-132
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Target Details	FABP4 encodes the fatty acid binding protein found in adipocytes. Fatty acid binding proteins are a family of small, highly conserved, cytoplasmic proteins that bind long-chain fatty acids and other hydrophobic ligands. It is thought that FABPs roles include fatty acid uptake, transport, and metabolism.
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.